




Patient Sling

SLING FOR ELOOR LIFTS

Seat/Back Sling




U-Sling w/Head Support



Sling w/Head Support & Insert Pocket




Seat/Back Sling w/Commode Opening




U-Sling w/o Head Support




Commode Sling w/Head Support & insert Pocket



Full Body Sling




One Piece Sling w/Positioning Strap



Sling w/Head Support




Full Body Sling w/Commode Opening



One Piece Commode Sling w/Positioning Strap



Commode Sling w/Head Support



Divided Leg Sling




Commode Seat and Back



SLING FOR STAND ASSIST LIFTS



Stand-Up Sling



Transfer Sling







Patient Sling











Patient Sling

No.	Item No.	Product Image	Specification
1	US001		<p>Face cloth: Polyester Filler: Resin Cotton Ribbon: Nylon Loading: 200kg Size: A*B*C mm S- 900*1180*720 M- 940*1260*780 L- 1040*1340*860</p> 
2	US002		<p>Face cloth: Polyester Filler: Resin Cotton Ribbon: Nylon Loading: 300kg Size: A*B*C mm S- 950*1220*790 M- 1020-1310*860 L- 1130*1430*950</p> 







Patient Sling

3	US003		<p>Face cloth: Polyester Filler: Sponge Ribbon: Nylon Loading: 220kg Size: A*B*C mm S- 830*1190*570 M- 900*1270*600 L- 1000*1380*680</p> 
4	US004		<p>Face cloth: Polyester Filler: EVA Ribbon: Nylon Loading: 200kg Size: A*B*C mm S- 820*930*860 M- 890*1000*920 L- 970*1070*1000</p> 
5	US005		<p>Face cloth: Polyester Filler: Sponge Ribbon: Nylon Loading: 200kg Size: A*B*C mm S- 870*930 M- 950*1020 L- 1040*1120</p> 








Patient Sling

6	US006		<p>Face cloth: Polyester Filler: EVA Ribbon: Nylon Loading: 220kg Size: A*B*C mm S- 830*1050*530 M- 910*1120*580 L- 1090*1400*720</p> 
7	US007		<p>Face cloth: Polyester Filler: Sponge Ribbon: Nylon Loading: 220kg Size: A*B*C mm S- 740*840*300 M- 830*920*330 L- 920*980*350</p> 



Patient Sling

8	US008		<p>Face cloth: Polyester Filler: Sponge Ribbon: Nylon Loading: 200kg Size: A*B*C mm S- 900*1020*300 M- 950*1100*320 L- 1130*1360*410</p> 
9	US009		<p>Face cloth: Polyester Filler: Sponge/Resin Cotton Ribbon: Nylon Loading: 225kg Size: A*B*C mm 1100*380*1230</p> 
10	US010		<p>Face cloth: Polyester Filler: Sponge /EVA Ribbon: Nylon Loading: 225kg Size: A*B*C mm 1000*320*1330</p> 