

Suction Cup Grab Bars

2013


FEATURE:

- * No drilling or screwing necessary, quickly installed on virtually or horizontally with no tools
- * Large release switches enable simple installation and removal.
- * Get a safe hold in the entire bathroom while taking a shower and ease to take while traveling.
- * CE Marked; SGS testing; FDA approval; Environmental friendly material-PAHs, Lead free.....etc.
- * Multi-national Patented safety indicators could timely detect the rail's adhesiveness




PATENTED SUCTION INDICATOR

"GREEN" means the correct attaching, otherwise, it will turn into the "RED" request to check the mounting surface and re-attach again. Do not place over joints or grout ines.







US
GERMANY
JAPAN
SPAIN
TAIWAN
CHINA







No.	Item No.	Product Image	Specification
1.	UNI-100 (NEW)		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 164 mm S- Cup Diameter: 93mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup <p>Packaging:</p> <ol style="list-style-type: none"> 1pc/inner box, 12 pcs/ctn Carton: 8PCS/Carton inner box : (L) 135 x (W) 135 x(H)145 mm carton box : (L) 430 x (W) 290 x (H) 315 mm

<p>2</p>	<p>UNI-100F (NEW)</p>		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 180 mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup <p>Packaging:</p> <ol style="list-style-type: none"> 1pc/inner box, 8pcs/ctn Carton: 8PCS/Carton inner box : (L) 145 x (W) 145 x(H)145 mm carton box : (L) 460 x (W) 290 x (H) 315 mm
<p>3</p>	<p>UNI-350</p>		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 400mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup Stainless Springs / screws <p>Packaging:</p> <ol style="list-style-type: none"> Inner box: 105 (L) x 110 (W) x 410 (H) mm Carton: 425 (L) x 235 (W) x 380 (H) mm 1pc/inner box, 6pcs/ctn N.W: 3.48 Kgs G.W: 4.82 Kgs
<p>4</p>	<p>UNI-400X</p>		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 436 mm S- Cup Diameter: 93 mm <p>Materials: Same as above</p> <p>Packaging:</p> <ol style="list-style-type: none"> Inner box:110(L) x 110(W) x 440(H) mm Carton: 345(L) x 225(W) x 455(H) mm 1pc/inner box, 6pcs/ctn N.W: 3.24 Kgs G.W: 4.55 Kgs
<p>5</p>	<p>UNI-500</p>		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 488mm S- Cup Diameter: 93 mm <p>Materials: Same as above</p> <p>Packaging:</p> <ol style="list-style-type: none"> Inner box: 97 (L) x 92 (W) x 495 (H) mm Carton: 507 (L) x 213 (W) x 300 (H) mm 1pc/inner box, 6pcs/ctn N.W: 3.24 Kgs G.W: 4.62 Kgs

6.	UNI-500X		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 545 mm S- Cup Diameter: 93 mm <p>Materials: Same as above</p> <p>Packaging:</p> <ol style="list-style-type: none"> Inner box:110(L) x 110(W) x 543(H) mm Carton: 345(L) x 225(W) x 560(H) mm 1pc/inner box, 6pcs/ctn N.W: 3.6 Kgs G.W: 5.36 Kgs
7	UNI-701 (NEW)		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 386mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup Stainless Springs / screws <p>Packaging:</p> <ol style="list-style-type: none"> Inner box:370(L) x 120(W) x 105(H) mm Carton: 385(L) x 260(W) x 335(H) mm 1pc/inner box, 6pcs/ctn
8	UNI-702 (NEW)		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 526mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup Stainless Springs / screws <p>Packaging:</p> <ol style="list-style-type: none"> Inner box:505(L) x 120(W) x 105(H) mm Carton: 520(L) x 260(W) x 335(H) mm 1pc/inner box, 6pcs/ctn
9	UNI-703 (NEW)		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 660mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup Stainless Springs / screws <p>Packaging:</p> <ol style="list-style-type: none"> Inner box:640(L) x 120(W) x 105(H) mm Carton: 655(L) x 260(W) x 335(H) mm 1pc/inner box, 6pcs/ctn

<p>10</p>	<p>UNI-704 (NEW)</p>		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 795mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup Stainless Springs / screws <p>Packaging:</p> <ol style="list-style-type: none"> Inner box:775(L) x 120(W) x 105(H) mm Carton: 790(L) x 260(W) x 335(H) mm 1pc/inner box, 6pcs/ctn
<p>11</p>	<p>UNI-705 (NEW)</p>		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 929mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup Stainless Springs / screws <p>Packaging:</p> <ol style="list-style-type: none"> Inner box:905(L) x 120(W) x 105(H) mm Carton: 920(L) x 260(W) x 335(H) mm 1pc/inner box, 6pcs/ctn
<p>12</p>	<p>UNI-801 (NEW)</p>		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 424mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup Stainless Springs / screws <p>Packaging:</p> <ol style="list-style-type: none"> Inner box: 410(L) x 120(W) x 105(H) mm Carton: 425(L) x 260(W) x 335(H) mm 1pc/inner box, 6pcs/ctn
<p>13</p>	<p>UNI-802 (NEW)</p>		<p>Specification:</p> <ol style="list-style-type: none"> Overall Length: 524mm S- Cup Diameter: 93 mm <p>Materials:</p> <ol style="list-style-type: none"> ABS hand rail PVC Suction Cup Stainless Springs / screws <p>Packaging:</p> <ol style="list-style-type: none"> Inner box: 510(L) x 120(W) x 105(H) mm Carton: 525(L) x 260(W) x 335(H) mm 1pc/inner box, 6pcs/ctn

14	UNI-803 (NEW)		Specification: 1. Overall Length: 734mm 2. S- Cup Diameter: 93 mm Materials: 1. ABS hand rail 2. PVC Suction Cup 3. Stainless Springs / screws Packaging: 1. Inner box: 715(L) x 120(W) x 105(H) mm 2. Carton: 730(L) x 260(W) x 335(H) mm 3. 1pc/inner box, 6pcs/ctn
15	UNI-804 (NEW)		Specification: 1. Overall Length: 834mm 2. S- Cup Diameter: 93 mm Materials: 1. ABS hand rail 2. PVC Suction Cup 3. Stainless Springs / screws Packaging: 1. Inner box: 815(L) x 120(W) x 105(H) mm 2. Carton: 830(L) x 260(W) x 335(H) mm 3. 1pc/inner box, 6pcs/ctn
16	UNI-805 (NEW)		Specification: 1. Overall Length: 934mm 2. S- Cup Diameter: 93mm Materials: 3. ABS hand rail 4. PVC Suction Cup 5. Stainless Springs / screws Packaging: 1. Inner box: 915(L) x 120(W) x 105(H) mm 2. Carton: 930(L) x 260(W) x 335(H) mm 3. 1pc/inner box, 6pcs/ctn


Obelis s.a.

Boulevard Général Wahis 53
1030 Brussels, Belgium
Tel : +(32) 2.732.59.54
Fax : +(32) 2.732.60.03
E-mail : mail@obelis.net